

MHC Staining Protocol

Procedure for Whole blood

- 1. Collect the blood into a tube containing anti-coagulant.
- 2. Lyse red blood cell using commercial red blood cell lysis reagent.
- 3. After red blood cell lysis, resuspend cells in FACS buffer (4% FBS in PBS) to a concentration of 1E7 cells/mL.
- 4. Add Dasatinib to the cells to a final concentration of 50nM and incubate 30 minutes.
- 5. Aliquot the cells into 100µL for each MHC sample. Add 1µg MHC tetramer and incubate for 30 minutes avoiding light.
- 6. If necessary, add any additional antibody (e.g anti-CD8) for cell sorting. Incubate for 30 minutes, avoiding light.
- 7. Centrifuge and remove the supernatant. Wash cells with FACS buffer.
- 8. Centrifuge and resuspend the cell pellet in 300µL FACS buffer and analyze by flow cytometry.

Procedure for Cell Preparations and Cell Suspensions

- 1. Prepare peripheral blood mononuclear cells (PBMC) or prepare a single-cell suspension from mouse spleen.
- 2. Centrifuge and resuspend the cell in FACS buffer (4% FBS in PBS). Adjust the cell concentration to 1E7 cells/mL in FACS buffer.
- 3. Add 50nM final concentration of Dasatinib into the cells and incubate for 30 minutes.
- 4. Aliquot the cells into 100µL for each MHC sample. Add 1µg MHC tetramer and incubate for 30 minutes, avoiding light.
- 5. If necessary, add any additional antibody (e.g anti-CD8) for cell sorting. Incubate for 30 minutes avoiding light.
- 6. Centrifuge and remove the supernatant. Wash the cell with FACS buffer.
- 7. Centrifuge and resuspend the cell pellet in 300µL FACS buffer and analyze by flow cytometry.

Additional Step: Peptide Pulsing

- 1. Prepare PBMC or mouse spleen cell under sterile conditions.
- 2. Using complete medium (RPMI 1640, 10% FBS, 500 IU/mL IL-2, 50 μ M β -Me, 1% P/S) dilute the cells to 2E6 cells/mL.
- 3. Dissolve the lyophilized peptide in DMSO at 10mg/mL. Add freshly dissolved peptide to a final concentration 10µg/mL in the medium. Incubate cells at 37°C for 72 hours.
- 4. After 72 hours incubation, feeding fresh complete medium every other day, cells can be grown for up to 2 weeks. The cells can be harvested for MHC staining after 7 days incubation, based on the cell counts.