




Introduction

Multi-pass transmembrane proteins (MP-TMPs) such as G protein-coupled receptors (GPCRs), ion channels, and transporters, play essential roles in maintaining normal physiological functions including signal transduction, cellular communication, and tissue homeostasis. Mutations in these proteins are frequently implicated in a wide range of diseases, making them highly attractive and validated targets for therapeutic antibody discovery.

However, producing stable and functionally active MP-TMPs remains a major bottleneck in drug research and development. Their reliance on a membrane environment for correct folding and activity poses significant challenges for recombinant expression and structural analysis. To address this unmet need, we leverage virus-like particles (VLPs) and nanodiscs as innovative platforms for the functional display of MP-TMPs in their native conformation. These systems provide native lipid bilayers that stabilize target proteins and enable downstream applications such as animal immunization, antibody screening, and analytical method development. Our approach offers a robust and scalable solution to accelerate antibody drug discovery against this critical class of targets.


Full-length STEAP1 proteins displayed on VLPs

Virus-like particles (VLPs) are non-infectious particles that mimic the structure of viruses but do not contain genetic material. They are often used as a tool of presenting multi-transmembrane proteins for various research purposes. KACTUS has successfully displayed STEAP1, an emerging new drug target for prostate cancer, as well as other types of multi-pass transmembrane proteins in a full-length, native conformation useful for boosting immune responses.



Full-length A2AR assembled into copolymer Nanodiscs

Nanodiscs have emerged as a powerful tool in functional and structural studies for membrane proteins. KACTUS nanodiscs are produced using SMA in a mammalian-cell-based, detergent-free process. The transmembrane segments are stabilized in the center of the phospholipid bilayer, surrounded by SMA, with intracellular and extracellular domains exposed. ARAR, a critical GPCR target for solid tumor indications, is soluble in nanodisc format, a native-like bilayer environment that maintains the physiological function of GPCRs and other transmembrane proteins.




References & Acknowledgements

- Pettersen, J. M., Yang, Y., & Robinson, A. S. (2023). Advances in nanodisc platforms for membrane protein purification. Trends in Biotechnology, 41(8), 1041-1054. <https://doi.org/10.1016/j.tibtech.2023.02.006>
- Tucker, D. F., Sullivan, J. T., Mattia, K., Fisher, C. R., Barnes, T., Mabila, M. N., Wilf, R., Sulli, C., Pitts, M., Payne, R. J., Hall, M., Deng, X., Davidson, E., Willis, S. H., Doranz, B. J., Chambers, R., & Rucker, J. B. (2018). Isolation of state-dependent monoclonal antibodies against the 12-transmembrane domain glucose transporter 4 using virus-like particles. Proceedings of the National Academy of Sciences, 115(22), E4990-E4999. <https://doi.org/10.1073/pnas.1716788115>

We would like to thank GatorBio (www.gatorbio.com) for their invaluable support of performing BLI assays on nanodiscs and sharing the data with KACTUS throughout their research.


VLP-displayed immunogens yielded functional antibodies with high affinity

KACTUS initiated an in-house immunization campaign using a CXCR4 Virus-like Particle (VLP) to evaluate its effectiveness in robust antiserum generation and to functionally validate the lead antibodies through ligand blocking assays. CXCR4 is a G protein-coupled receptor (GPCR) that plays a significant role in cell signaling, leukocyte migration, and development. CXCR4 is involved in a series of oncological and autoimmune diseases, making it a prominent target for antibody development.




Biotinylated VLPs & Nanodiscs for phage panning

Phage panning usually requires biotinylated proteins for optimal binding and screening. To fulfill this need, we developed biotinylated VLPs and Nanodiscs with full-length proteins displayed on their surface without interfering with the protein conformation. We showed that both our biotinylated VLP and Nanodiscs produced robust ELISA signals, demonstrating their ability to support phage panning, a screening process requiring strong binding interactions for effective target enrichment.



Biologically-active TCR-CD3 complex Nanodisc for CD3 antibody screening

Traditionally, anti-CD3 antibodies are generated and screened using soluble CD3 subunits. These often show limited effectiveness as CD3 proteins naturally exist as part of the full TCR-CD3 complex on the T cell surface. Antibodies raised against individual subunits may fail to bind the native structure on the cell surface. To improve biological relevance and functionality, KACTUS has successfully produced nanodiscs containing the full complement of TCR-CD3 subunits, including 2XCD3ε, CD3γ, CD3δ, 2XCD3ζ, TCRα, and TCRβ, to facilitate antibody screening and lead isolation.



Conclusion

- KACTUS VLP and SMA-based nanodisc display platform allows robust, soluble display of membrane proteins with native structural conformation and robust bioactivity.
- KACTUS VLP and Nanodisc-displayed antigens are powerful tools for antibody discovery and bioanalytical assay development against membrane proteins, which can be applied for animal immunization, ELISA, SPR and BLI assays.

Membrane Proteins produced by KACTUS

We have successfully produced a variety of membrane proteins including GPCRs, Tetraspanins and GIPRs. To request a sample or learn more about our custom production capability, please contact us at support@kactusbio.us.