




Custom Protein Production Services



Tailored production of complex & specialized proteins


Complete expertise in construct design, expression characterization, and structure optimization



6-8 Weeks Turnaround Time



Proteins from Multiple Species



Service Features

Example Proteins

MHC Tetramers
Soluble TCRs
Transmembrane Proteins
- GPCRs
- Solute Carriers
- 4HB/Tetraspanins
- Ion Channels
Other Enzymes & Proteins

Expression Systems


Mammalian
Bacterial
Yeast
Insect

Modifications


Tag: His, hFc, Avi, etc.
Species: Human, Cyno, Mouse, etc.
Biotinylation
Fluorescent labeling
Full-length protein or single domain
Liquid or lyophilized



Expression
Optimization



Custom Purity and
Endotoxin Levels



Bioactivity
Testing



Flexible
Production Scales



Transparent
Collaboration

Learn more
& request quote:

